

HEMO Weld



DELCON
Our technology. Your health.




WE.

KNOW HOW

HEMOWELDs are designed to produce sterile seals on standard PVC tubing normally used in conjunction with blood bags.

You can download full specifications and the Technical Data Sheet from our website www.delcon.it



HEMO Weld-B

Portable tube sealer battery powered

- ✓ Lightweight, about 1 kg
- ✓ Powerful, 500 sealing procedures
- ✓ Certified to work with needle still connected to the patient
- ✓ Supplied with a shoulder bag

HEMOWELD-B is supplied together with a comfortable shoulder bag. It is portable and equipped with a powerful rechargeable battery pack.

It can be used to seal most of standard tubing available on the market. It is easy to use and to clean and thanks to its low power consumption it permits more than 500 procedures with fully charged battery. Each seal appears with a central line to facilitate the snap-apart separation of the two parts.

An anti-spark and overheat protection systems to avoid damage to the system and tubing it is included as standard. A beep and LED signals indicate the presence of sparks while sealing, dirty electrodes or low battery. The sealer is certified as CF Class for the Electrical Safety Standards and it can be used with patient/donor still connected to the needle.

ORDERING INFORMATION

HWELD-B.DB.00

HEMOWELD-B kit including RF Generator, battery charger, shoulder bag, spare coaxial cable

OPTIONAL ACCESSORY

HWELDB.DB.30

External battery charger 100-240V for Hemoweld-B



HEMO Weld-T

Automatic bench top sealer for PVC tubing

- ✓ For heavy duty job
- ✓ Automatic clamp
- ✓ Self adjusting RF energy
- ✓ Optional handle for best flexibility of use

This sealer is designed for bench use and for heavy duty job. The sealing clamp automatically closes once the tube is inserted and the seal is done in about 1.5 seconds. The protection cap of the clamp is easy to remove and it is washable. HEMOWELD-T can be used to seal most of standard PVC tubing available on the market.

It is easy to use and to clean and it has a high working capacity. Each seal appears with a central line to facilitate the snap-apart separation of the two parts: sterility of internal blood components is maintained as well. An anti-spark and overheat protection systems to avoid damage to the system and tubing it is included as standard.

A beep and LED signals indicate the presence of sparks while sealing, dirty electrodes or system over heating. A separate handle is available as an option to improve the flexibility of work.

ORDERING INFORMATION

HWELD-T.DB.00

Bench tube sealer

OPTIONAL ACCESSORY

SEAL-HAND1

Standard sealing handle equipped with two coaxial cables



HEMO Weld-MC

Automatic bench sealer for segments preparation

- ✓ System with four clamps
- ✓ Each clamp can be enabled or disabled
- ✓ Four segments prepared in few seconds
- ✓ Self adjusting RF energy
- ✓ Optional handle for best flexibility of use

HEMOWELD- MC is designed to rapidly prepare up to four segments of PVC tubing. Each segment has a standard length of 75 mm, enough to include any printed code with sufficient volume to perform any testing. The tubing is loaded inside a moving carousel and once correctly in place the sealing process starts automatically in sequence. The whole cycle with all four clamps enabled lasts about 7 seconds. The Operator can chose to use all four clamps or only some at his choice by a programming key. The protection cover of the clamps is easy to remove and it is washable.

HEMOWELD-MC can be used to seal most of standard PVC tubing available on the market. The sealer is easy to use and to clean and it has a high working capacity. Each seal appears with a central line to facilitate the snap-apart separation of the two parts: sterility of internal blood components is maintained as well.

An anti-spark and overheat protection systems to avoid damage to the system and tubing it is included as standard. A beep and LED signals indicate the presence of sparks while sealing, electrodes dirty or system over heating. A separate handgun is available as an option to improve the flexibility of work.

ORDERING INFORMATION

HWELD-MC.DB.00

Bench sealer for segments preparation

OPTIONAL ACCESSORY

SEAL-HAND1

Standard sealing handle equipped with two coaxial cables





DELCON

Our technology Your health

DELCON srl

20862 Arcore (MB) - Italy - Via Matteucci, 25/27

Phone +39 039.617670

Fax +39 039.617677

www.delcon.it



Delcon reserves the right to modify, add, and remove the functional and visual features of the devices, at its sole discretion, at any time, and without prior notice

M A D E I N I T A L Y

CE



Management System
ISO 9001:2008
EN ISO 13485:2012
www.tuv.com
ID: 9106018273

